

Vascular Intervention // Coronary Drug-Eluting Stent System



Orsiro® Mission des

Even better deliverability for the outstanding Orsiro DES



Orsiro Mission DES Even better deliverability for the outstanding Orsiro DES

The next level of deliverability¹

1st in Push⁴

Transmitting up to 57% more force from hub to tip.



1st in Track⁴

Up to **30% less** force needed to follow the path to the lesion.



1st in Cross⁴

Up to 75% less force needed to successfully cross demanding anatomies.







"Lesion crossing with low friction, reliable performance"

Dr. Mathias Brandt, Paracelsus Medical University, Salzburg, Austria



3 out of 4

cases show Orsiro Mission DES to have better deliverability in complex lesions $\Delta 5$



* ø 2.25 – 3.0 mm

△ 'Complex lesions' are defined by lesions where both calcification and tortuosity

were individually reported as 'Moderate or Severe'. In comparison to Xience Sierra and Synergy.



with kink resistance

 \bigcirc

Ultrathin struts²

Strut thickness in perspective⁷

For early endothelialization

Orsiro BIOTRONIK CoCr-SES OCr-SES OCr-SES OCr-SES

60 µm*

Synergy Boston Scientific PtCr-EES

74 µm



Strut coverage¹⁰

90 days[∆]









\bigcirc

Outstanding patient outcomes³

One of the most studied DES[§]

More than 55,000 patients enrolled and more than 68 Studies started^{ϕ} BIOFLOW-V, FDA pivotal trial (n = 1,334)^{12, 13, 14, 15}

Orsiro Mission DES is indicated for complex patients and lesions, including:*

BIOSTEMI (n=1,300)

Continued Superiority in STEMI at 2 years.^{¶16}

Target Lesion Failure (TLF) rate at 2 years.

Rate Ratio (95% BCI**): 0.58 (0.40-0.84) Posterior probability of Superiority: 99.8% Bayesian ITT Population[¤]

BIO-RESORT Small Vessels (n=1,506)

Target Lesion Revascularization (TLR) rate at 3 yrs.¹⁷

- § in large RCTs, based on Taglieri et al. Meta-analysis, against currently used DES
 \$ BIOTRONIK data on file, status January 2020
- [¶] vs. Xience, based on TLF, in the BIOSTEMI trial
- ◊ p-values for 36-m frequentist analysis (see supplemental material).
- *As per IFU: ACS Acute Coronary Syndrome; STEMI ST-Elevation Myocardial Infarction; DM – Diabetes Mellitus.
- HBR High Bleeding Risk; B2C Complex Lesions; SV Small Vessels; MVD Multi-Vessel Disease. ** BCI: Bayesian Credibility Interval.
- ¤ n= 1,300 newly enrolled STEMI patients including 407 patients from the BIOSCIENCE STEMI subgroup used as prior information.
- [‡] Based on a Rate Ratio of 0.58.

Orsiro® Mission des

The Orsiro Mission Sirolimus-Eluting Coronary Stent System is a drug-eluting balloon-expandable stent pre-mounted on a rapid-exchange PTCA catheter delivery system. Vascular Intervention Coronary

Indication				
	Orsiro Mission DES is indicated for improving coronary luminal diameter in patients with symptomatic ischemic heart disease due to discrete de-novo stenotic lesions and in-stent restenotic lesions (length ≤ 40 mm) in the native coronary arteries with a reference vessel diameter of 2.25 mm to 4.0 mm including the following patient and lesion subsets:			
	Acute Coronary Syndrome (ACS) ST-Elevation Myocardial Infarction Diabetes Mellitus (DM) Complex Lesions (B2/C) High Bleeding Risk (HBR)	Long Lesions (LL) (e.g. ≥ 20 mm) (STEMI) Small Vessels (SV) (e.g. ≤ 2.75 mm) Multi-Vessel Disease (MVD) Male/Female Old Patients (e.g. > 65 y)		
Technical Data	Stent			
	Stent material (Cobalt chromium, L-605		
	Strut thickness	ø 2.25 – 3.0 mm: 60 μm (0.0024"); ø 3.50 – 4.0 mm: 80 μm (0.0031")		
	Passive coating	roBIO (Amorphous Silicon Carbide)		

<u> </u>			
Active coating	BIOlute bioabsorbable Poly-L-Lactide (PLLA) eluting		
	a limus drug		
Drug dose	1.4 µg/mm²		

Delivery system

Catheter type	Rapid exchange
Recommended guide catheter	5F (min. I.D. 0.056")
Guide wire diameter	0.014"
Usable catheter length	140 cm
Balloon material	Semi crystalline polymer
Coating (Distal shaft)	Hydrophilic
Coating (Proximal shaft)	Hydrophobic
Marker bands	Two swaged platinum-iridium markers
Lesion entry profile	0.017"
Distal shaft diameter	2.7F: ø 2.25 – 3.0 mm; 2.9F: ø 3.5 – 4.0 mm
Proximal shaft diameter	2.0F
Nominal pressure (NP)	10 atm
Rated burst pressure (RBP)	16 atm

Storage

Use Before Date (UBD)	24 months
Temperature	Between 15°C (59°F) and 25°C (77°F), short term excursions between 10°C (50°F) and 40°C (104°F) are allowed

35	40
419143	419149
419144	419150
419145	419151
419146	419152
419147	419153
419148	419154
	35419143419144419145419146419147419148

1. In comparison to Xience Sierra, Resolute Onyx and Synergy for bench tests on pushability, trackability and crossability, BIOTRONIK data on file; 2. As characterized with respect to strut thickness in Bangalore et al. Meta-analysis; 3. Based on investigator's interpretation of BIOFLOW-V primary endpoint result; 4. BIOTRONIK data on file; 5. Evaluation of Market Acceptance, BIOTRONIK data on file; 6. Per investigators' interpretation of preclinical studies with Orsiro as mentioned in Cassese et al. J Thorac Dis 2018;10(2):688-692; 7. Stefanini GG et al. Coronary stents: novel developments. Heart. 2014 Jul 1;100(13):1051-61; 8. Low AF. Stent platform for procedural success: Introducing the Continuous Sinusoidal & Core Wire Technologies. Presented at: AsiaPCR; 22-24 January, 2015; Singapore, Singapore; 9. Tolentino A. Evolving DES Strategy: Biodegradable Polymervs. Bioabsorbable Scaffold. Presented at: Cardiovascular Nurse/Technologist Symposium; June 17, 2016; New York, USA; 10. Secco G et al. Time-related changes in neointimal tissue coverage of a novel Sirolimus eluting stent: Serial observations with optical coherence tomography. Cardiovascular Revascularization Medicine 17.1 (2016): 38-43; 11. Pilgrim T et al. 5-year outcomes of the BIOSCIENCE randomised trial. Supplementary appendix; Lancet 2018; published online Aug 28. http://dx.doi.org/10.1016/ S0140-6736(18)31715-X; 12. Kandzari D, et al. BIOFLOW-V: A Prospective Randomized Multicenter Study to Assess the Safety and Effectiveness of the Orsiro Sirolimus Eluting Coronary Stent System in the Treatment Of Subjects With up to Three De Novo or Restenotic Coronary Artery Lesions Science. Presentation at E SC 2017; 13. Kandzari D et al. Ultrathin Bioresorbable Polymer Sirolimus-Eluting Stents versus Thin Durable Polymer Everolimus-Eluting Stents: Journal of American College of Cardiology (2018), doi: https://doi.org/10.1 016/j.jac c.2018. 09.019; 14. Kandzari D et al. J Am Coll Cardiol. Cardiovasc Interven. 2020, doi: 10.1016/ j.jcin.2020.02.019; 15. Kandzari D et al. J Am Coll Cardiol. Cardiovasc Interven. 2020. Supplemental Material; 16. Pilgrim et al. Biodegradable - versus durable-polymer drug-eluting stents for STEMI. Final 2-year outcomes of the BIOSTEMI trial. J Am Coll Cardiol. Cardiovasc Interven. 2021, doi: 10.1016/j.jcin.2020.12.011; 17. Buiten R et al. Outcomes in patients treated with thin-strut, very thin-strut, or ultrathin-strut drug-eluting stents in small coronary vessels – A prespecified analysis of the randomized BIO-RESORT trial; JAMA Cardiol. Published online May 21, 2019. doi:10.1001/jamacardio.2019.1776; ClinicalTrials.gov: NCT01674803. Orsiro, Orsiro Mission, proBIO and BIOlute are trademarks or registered trademarks of the BIOTRONIK Group of Companies. Synergy and Promus are trademarks or registered trademarks of the Boston Scientific Group of Companies. Resolute, Resolute Onyx and Integrity are trademarks or registered trademarks of the Medtronic Group of Companies. Xience and Xience Sierra are trademarks or registered trademarks of the Abbott Group of Companies. Ultimaster is a trademark or registered trademark of the Terumo Group of Companies. BioMatrix is a trademark or registered trademark of the Biosensors International Group.

Clinical data conducted with Orsiro, Orsiro Mission's predecessor device can be used to illustrate Orsiro Mission clinical outcomes.

BIOTRONIK AG Ackerstrasse 6 8180 Bülach, Switzerland Tel +41 (0) 44 8645111 Fax +41 (0) 44 8645005 info.vi@biotronik.com www.biotronik.com

© 2021 BIOTRONIK AG – All rights reserved. Specifications are subject to modification, revision and improvement.

